

# ELISA Kit for Heat Shock 70kDa Protein 1A (HSPA1A) Catalog No: tcue148

**Available Sizes** 

**Size:** 96T

Specifications

## **Research Area:**

Signal transduction; Tumor immunity; Infection immunity; Developmental science;

**Species Reactivity:** Sus scrofa; Porcine (Pig)

# Sample Type:

Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids

## Sensitivity:

The minimum detectable dose of this kit is typically less than 0.61ng/mL

#### **Detection Range:**

1.56-100ng/mL

# Assay Time:

3h

## **Detection Method:**

Enzyme-linked immunosorbent assay for Antigen Detection.

# **Tested Application:**

ELISA

# **Test Principle**

The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Heat Shock 70kDa Protein 1A (HSPA1A). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Heat Shock 70kDa Protein 1A (HSPA1A). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Heat Shock 70kDa Protein 1A (HSPA1A), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a



change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm  $\pm$  10nm. The concentration of Heat Shock 70kDa Protein 1A (HSPA1A) in the samples is then determined by comparing the O.D. of the samples to the standard curve.



All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!

Copyright 2021 Taiclone Biotech Corp.