

# Jurkat, Clone E6-1

Catalog No: tcel129



## Available Sizes

**Size:** 1×10<sup>6</sup>cells/t25culturebottle



## Specifications

### Application:



transfection

### Subculturing:





Cultures can be maintained by the addition of fresh medium or replacement of medium. Alternatively, cultures can be established by centrifugation with subsequent resuspension at the subcultivation ratio of 1:2-1:4.

## Product Description

Clone E6-1 cells produce large amounts of IL-2 after stimulation with phorbol esters and either lectins or monoclonal antibodies against the T3 antigen(both types of stimulants are needed to induce IL-2 production.







All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!