

Jurkat, Clone E6-1

Catalog No: tcel129



Available Sizes

Size: 1×10⁶cells/t25culturebottle



Specifications

Application:



transfection





Subculturing:

Cultures can be maintained by the addition of fresh medium or replacement of medium. Alternatively, cultures can be established by centrifugation with subsequent resuspension at the subcultivation ratio of 1:2-1:4.

Product Description

Clone E6-1 cells produce large amounts of IL-2 after stimulation with phorbol esters and either lectins or monoclonal antibodies against the T3 antigen(both types of stimulants are needed to induce IL-2 production).





All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!