

## KG-1 Catalog No: tcel132



**Available Sizes** 

Size: 1×10<sup>6</sup> cells/t25 culture bottle

Specifications

## **Application:**

transfection

## Subculturing:

Cultures can be maintained by the addition of fresh medium or replacement of medium. Alternatively, cultures can be established by centrifugation with subsequent resuspension at the subcultivation ratio of 1:2-1:3.

## **Product Description**

KG-1 myeloblasts can differentiate into macrophages in the presence of phorbol esters. They show a good response to colony stimulating factor(CSF). The line is EBNA negative(EBNA-).





Copyright 2021 Taiclone Biotech Corp.



Web: www.taiclone.com Tel: +886-2-2735-9682 Email: order@taiclone.com



All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!

Copyright 2021 Taiclone Biotech Corp.