

## PC-12

**Catalog No: tcel480**



### Available Sizes

---

**Size:** 1×10<sup>6</sup>cells/t25culturebottle

---



### Specifications

#### **Subculturing:**

Remove and discard culture medium. Briefly rinse the cell layer with DPBS solution to remove all traces of serum that contains trypsin inhibitor. Add 1.0 to 2.0 mL of Trypsin-EDTA solution to flask and observe cells under an inverted microscope until cell layer is dispersed (usually within 2 to 3 minutes). Cells that are difficult to detach may be placed at 37°C to facilitate dispersal. Add 4.0 to 6.0 mL of complete growth medium and aspirate cells by gently pipetting. Add appropriate aliquots of the cell suspension to new culture vessels.

---

### Product Description

The cells do not synthesize epinephrine. This adherent variant(PC-12 Adh) has been adapted to Corning CellBIND<sup>®</sup> flasks to improve cell attachment.



All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!