

Elastic Stain Kit (Modified Verhoff's)

Catalog No: tcst1983



Available Sizes

Size: 100 tests



Specifications

Research Area:

Histochemistry

Component:

Hematoxylin Solution (5%) 250ml; Ferric Chloride (10%, Aqueous) 125ml ;Lugol's Iodine Solution 125ml ;Ferric Chloride (2%) Differentiating Solution 125ml ;Sodium Thiosulfate Solution (5%) 125ml ;Van Gieson's Solution 125ml

Storage Instruction:


18~25 °C Keep away from fire or flame.

Notes

Interpretation data: Elastic fibers: Black to Blue/Black Nuclei: Blue to Black Collagen: Red Muscle & Other: Yellow

Product Description

Taiclone Elastic Stain Kit is used for staining procedure histological of elastin in tissue sections. Demonstration of elastic tissue is useful in cases of emphysema (atrophy of elastic tissue), arteriosclerosis (thinning and loss of elastic fibers) and various other vascular diseases.



All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!