

ELISA Kit for Pregnancy Associated Plasma Protein A (PAPPA)

Catalog No: tcue81

Available Sizes

Size: 96T



Research Area: Reproductive science;Genetic science;

Species Reactivity:

Mus musculus (Mouse)

Sample Type: serum, plasma, tissue homogenates and other biological fluids

Sensitivity:

The minimum detectable dose of this kit is typically less than 29pg/mL

Detection Range:

78-5,000pg/mL

Assay Time: 3h

Detection Method:

Enzyme-linked immunosorbent assay for Antigen Detection.

Tested Application:

ELISA


SwissProt: Q8R4K8

Copyright 2021 Taiclone Biotech Corp.



Test Principle

The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Pregnancy Associated Plasma Protein A (PAPPA). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Pregnancy Associated Plasma Protein A (PAPPA). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Pregnancy Associated Plasma Protein A (PAPPA), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm \pm 10nm. The concentration of Pregnancy Associated Plasma Protein is then determined by comparing the O.D. of the samples to the standard curve.



All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!

Copyright 2021 Taiclone Biotech Corp.