



ELISA Kit for Neuronal Apoptosis Inhibitory Protein (NAIP)

Catalog No: tcue726


Available Sizes Size: 96T	
Research Area: Neuro science;	
Species Reactivi Homo sapiens (Hu	
Sample Type: Tissue homogenat	es and other biological fluids.
Sensitivity: The minimum dete	ectable dose of this kit is typically less than 0.058ng/mL
Detection Range 0.156-10ng/mL	}:
Assay Time: 3h	
Detection Metho Enzyme-linked im	od: munosorbent assay for Antigen Detection.
Tested Applicati ELISA	on:
SwissProt: Q13075	





Test Principle

The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Neuronal Apoptosis Inhibitory Protein (NAIP). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Neuronal Apoptosis Inhibitory Protein (NAIP). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Neuronal Apoptosis Inhibitory Protein (NAIP), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Neuronal Apoptosis Inhibitory Protein (NAIP) in the samples is then determined by comparing the O.D. of the samples to the standard curve.



All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!