



ELISA Kit for Glycoprotein 2, Zymogen Granule Membrane (GP2)

Catalog No: tcue3313


A	vailable Sizes
Size: 96T	
S _I	pecifications
Researc Endocrine	
	Reactivity: orvegicus (Rat)
Sample Serum, p	Type: lasma, tissue homogenates and other biological fluids.
Sensitiv The mini	ity: mum detectable dose of this kit is typically less than 0.122ng/mL
Detectio 0.312-20	on Range: ng/mL
Assay Ti 3h	ime:
	on Method: linked immunosorbent assay for Antigen Detection.
Tested A	Application:
SwissPro	ot:





Test Principle

The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Glycoprotein 2, Zymogen Granule Membrane (GP2). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Glycoprotein 2, Zymogen Granule Membrane (GP2). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Glycoprotein 2, Zymogen Granule Membrane (GP2), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Glycoprotein 2, Zymogen Granule Membrane (GP2) in the samples is then determined by comparing the O.D. of the samples to the standard curve.



All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!