

M-NFS-60 [NFS-60]

Catalog No: tcel301



Available Sizes

Size: 1×10⁶cells/t25culturebottle





Specifications





Subculturing:

Cultures can be maintained by the addition of fresh medium or replacement of medium. Alternatively, cultures can be established by centrifugation with subsequent resuspension at 1 to 2×10⁵ viable cells/ml. Maintain cell density between 1×10⁵ and 1×10⁶ viable cells/ml.

Product Description

The cells are responsive to both interleukin 3(interleukin-3, IL-3) and macrophage colony stimulating factor(M-CSF).The cells contain a truncated c-myb proto-oncogene caused by integration of a retrovirus.





All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!