

## M-NFS-60 [NFS-60]

Catalog No: tcel301



### Available Sizes

---

**Size:** 1×10<sup>6</sup>cells/t25culturebottle

---



### Specifications



#### Subculturing:

Cultures can be maintained by the addition of fresh medium or replacement of medium. Alternatively, cultures can be established by centrifugation with subsequent resuspension at 1 to 2×10<sup>5</sup> viable cells/ml. Maintain cell density between 1×10<sup>5</sup> and 1×10<sup>6</sup> viable cells/ml.





---

### Product Description

The cells are responsive to both interleukin 3(interleukin-3, IL-3) and macrophage colony stimulating factor(M-CSF).The cells contain a truncated c-myb proto-oncogene caused by integration of a retrovirus.







All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!