

## PLC / PRF / 5

### Catalog No: tcel415



#### Available Sizes

---

**Size:** 1×10<sup>6</sup>cells/t25culturebottle

---



#### Specifications

##### **Subculturing:**

Remove and discard culture medium. Briefly rinse the cell layer with DPBS solution to remove all traces of serum that contains trypsin inhibitor. Add 1.0 to 2.0 mL of Trypsin-EDTA solution to flask and observe cells under an inverted microscope until cell layer is dispersed (usually within 2 to 3 minutes). Cells that are difficult to detach may be placed at 37°C to facilitate dispersal. Add 4.0 to 6.0 mL of complete growth medium and aspirate cells by gently pipetting. Add appropriate aliquots of the cell suspension to new culture vessels.

---

#### **Product Description**

The line was originally contaminated with mycoplasma, and was cured by treatment with BM-cycline. The cells secrete HBsAg.



All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!