

RT4

Catalog No: tcel431



Available Sizes

Size: 1×10⁶cells/t25culturebottle



Specifications

Subculturing:

Remove and discard culture medium. Briefly rinse the cell layer with DPBS solution to remove all traces of serum that contains trypsin inhibitor. Add 1.0 to 2.0 mL of Trypsin-EDTA solution to flask and observe cells under an inverted microscope until cell layer is dispersed (usually within 2 to 3 minutes). Cells that are difficult to detach may be placed at 37°C to facilitate dispersal. Add 4.0 to 6.0 mL of complete growth medium and aspirate cells by gently pipetting. Add appropriate aliquots of the cell suspension to new culture vessels.

Product Description

A contamination with Mycoplasma orale was cured in October 1971.



All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!