



# L5178Y TK+ / - clone (3.7.2C)

**Catalog No: tcel455** 



### **Available Sizes**

**Size:** 1×10<sup>6</sup>cells/t25culturebottle



## **Specifications**

#### **Subculturing:**

Cultures can be maintained by the addition of fresh medium or replacement of medium. Alternatively, cultures can be established by centrifugation with subsequent resuspension at 1 to  $2\times10^5$  viable cells/ml. Maintain cell density between  $1\times10^5$  and  $1\times10^6$  viable cells/ml.

## **Product Description**

This line is a derivative of L5178Y-S



All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!