

## L5178Y TK+ / - clone (3.7.2C)

Catalog No: tcel455



### Available Sizes

---

**Size:** 1×10<sup>6</sup>cells/t25culturebottle

---



### Specifications

#### **Subculturing:**

Cultures can be maintained by the addition of fresh medium or replacement of medium. Alternatively, cultures can be established by centrifugation with subsequent resuspension at 1 to 2×10<sup>5</sup> viable cells/ml. Maintain cell density between 1×10<sup>5</sup> and 1×10<sup>6</sup> viable cells/ml.

---

### Product Description

This line is a derivative of L5178Y-S



All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!