

## **MNNG / HOS CI #5 [R-1059-D]**

**Catalog No: tcel492**



### **Available Sizes**

---

**Size:** 1×10<sup>6</sup>cells/t25culturebottle

---



### **Specifications**

#### **Subculturing:**

Remove and discard culture medium. Briefly rinse the cell layer with DPBS solution to remove all traces of serum that contains trypsin inhibitor. Add 1.0 to 2.0 mL of Trypsin-EDTA solution to flask and observe cells under an inverted microscope until cell layer is dispersed (usually within 2 to 3 minutes). Cells that are difficult to detach may be placed at 37°C to facilitate dispersal. Add 4.0 to 6.0 mL of complete growth medium and aspirate cells by gently pipetting. Add appropriate aliquots of the cell suspension to new culture vessels.

---

### **Product Description**

The cells exhibit a high saturation density, a high plating efficiency in soft agar and are tumorigenic in nude mice.



All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!