

BV2 Cells Complete Medium

Catalog No: tcem2474

Available Sizes

Size: 4x125ml



Component: MEM[10% FBS[0.1mg/mL Streptomycin+100U/mL Penicillin G

Endotoxin Level:

≤5

Storage Instruction:

2-8°C, Shading light for 3 months

Notes

1. This product is sterilized by 0.22 μ m filtration. 2. It is necessary to pay attention to the aseptic operation and avoid the pollution during the culture. 3. This product contains serum, penicillin and streptomycin. If there is no special need, don't add serum, penicillin and streptomycin. It can be used directly. 4. It is not suitable for long time storage at room temperature or 4 °C.

Product Description

Ready to use Media for BV2 Cells.



Copyright 2021 Taiclone Biotech Corp.



All products are for RESEARCH USE ONLY. Not for diagnostic & therapeutic purposes!

Copyright 2021 Taiclone Biotech Corp.